

M9 Minimal Salts 5x, Powder

Cat. No. 48505

Product Description:

General M9 Minimal Salts 5x are suitable for the preparation of M9 Minimal Medium for cultivation of recombinant *E. coli* strains.

Range of application of M9 Minimal Medium

- Maintaining positive selection pressure on plasmids coding for genes to produce essential substances such as amino acids and vitamins
- Maintenance of F⁺-containing bacteria strains for M13 propagation
- After supplementation with specific amino acids or other metabolites for selection of specific auxotrophs

Composition M9 Minimal Salts 5x is a fivefold concentrate.

Component	Concentration
Na ₂ HPO ₄	10 g/l
KH ₂ PO ₄	5 g/l
NH ₄ Cl	5 g/l
NaCl	2.5 g/l

Storage Recommended storage temperature of powder is +15 °C – +30 °C. Keep container tightly closed, because powder is very hygroscopic.

Solution procedure:

1. For 1 liter 5x concentrate solve 52.2 g powder in 1 l dest. water.
2. For sterilisation autoclave for 15 min. at 121 °C.

Preparation of M9 minimal medium (example, composition application-specific):

1. Add 200 ml sterile 5x M9 Minimal Salt solution to 750 ml sterile dest. H₂O (cooled to 45 – 50 °C), adjusting the final volume to 1 liter.
 2. Aseptically add 20 ml filter-sterilized 20 % glucose solution, 2 ml sterile 1 M MgSO₄ solution and, if desired, 0.1 ml sterile 1 M CaCl₂ solution. Mix well.
 3. If desired, supplement with amino acids, as appropriate.
-

¹ Sambrook, J. et al., Molecular Cloning : A Laboratory Manual, 2nd ed., p. A.3, Cold Spring Harbor Laboratory Press, Cold Spring Harbor, NY.

Version 03/07