

# Mammalian Protein Extraction Kits

**SERVA**  
■ serving scientists ■

## Fast and Simple Isolation of Native Nuclear, Cytoplasmic and Membrane Proteins

The mammalian protein extraction kits provide a fast and efficient method to extract total proteins or cytoplasmic, nuclear and membrane protein fractions from mammalian cells or tissues.

The lysis buffers allow for mild but efficient lysis with high protein yield. The performance is optimized and validated for several cell lines and tissues. There is no need for mechanical cell disruption. The single formulation reagents provide versatility, because protease or phosphatase inhibitors, reducing agents etc. may be added prior to cell lysis. The reagents are compatible with standard protein quantification assays like BCA Protein Assay.



- Contain ready-to-use optimized extraction buffers and protease inhibitor cocktail
- High protein yields from cells or tissues
- Gentle formulation preserves protein activity for downstream assays
- No need for mechanical disruption or ultracentrifugation
- Lysates are directly compatible with immunoprecipitation, SDS-PAGE, Western Blot, ELISA, EMSA and other functional assays

### Low Cross Contamination




Fig. 1 Membrane proteins were isolated from HEK293 cells using the SERVA Mammalian Membrane Protein Extraction Kit and the Membrane Protein Extraction Kit of a competitor, following the respective protocols. To determine the percentage of cross contamination by cytoplasmic proteins in the membrane protein fraction,  $\beta$ -Tubulin was detected in the cytoplasmic (lane 1,2) and membrane extract (lane 3,4) by Western blot.




Fig. 2 Cytoplasmic and nuclear proteins were isolated from HEK293 cells using the SERVA Mammalian Membrane Protein Extraction Kit and the Membrane Protein Extraction Kit of a competitor, following the respective protocols. To analyze the efficiency of cellular fractionation,  $\beta$ -Tubulin and HDAC1 were detected in the cytoplasmic (lane 1,2) and nuclear extract (lane 3,4) by Western blot.

### High Protein Yield from Cells and Tissues




Fig. 3 Total proteins were isolated from HeLa cells and A549 cells and rat lung tissue ( $1 \times 10^7$  cells and 50 mg tissue) using the SERVA Mammalian Total Protein Extraction Kit. 20  $\mu$ g of the extract were run on a 4-20 % SDS PAGE gel and analyzed by Western Blot using 1:200 dilutions of primary antibodies to various specific proteins.

### Extracted Proteins Are Directly Compatible with Enzyme Activity Assays




Fig. 4 Enzymatic activity assay of Na<sup>+</sup>/K<sup>+</sup>-ATPase  $\alpha$  (lane 2, lane 1 = negative control) shows that proteins isolated with the Mammalian Protein Extraction Kit are fully functional and that activity is preserved for downstream assays.

# Mammalian Protein Extraction Kits


**SERVA**  
■ serving scientists ■

## Mammalian Total Protein Extraction Kit




Extract total native proteins (cytoplasmic, membrane and nuclear proteins) from mammalian cells and tissues.

## Mammalian Membrane Protein Extraction Kit



Extract native membrane and cytoplasmic proteins from mammalian cells and tissues. Efficiency is up to 90 % for membrane proteins with at least 1 – 2 trans-membrane domains. Cross-contamination of cytoplasmic proteins is usually less than 10 %.

## Mammalian Nuclear and Cytoplasmic Protein Extraction Kit



Extract native nuclear and cytoplasmic proteins from mammalian cells and tissues. Cross-contamination of cytoplasmic proteins is usually less than 10 %.

## Ordering Information

Product	Size	Cat. no.
Mammalian Total Protein Extraction Kit	100 Samples	39241.01
Mammalian Nuclear and Cytoplasmic Protein Extraction Kit	50 Samples	39243.01
Mammalian Membrane Protein Extraction Kit	50 Samples	39242.01