

Albumin bovine modified Cohn Fraction V, pH 5.2
standard grade, lyophil.

Cat.No. : 11941

Parameter	Method	Specification
Molecular weight		ca. 67 000
Appearance	Visual inspection	white to yellow with tan to green cast lyophilisate
Solubility	4% in H ₂ O	clear to slightly hazy, faint yellow to yellow-green solution
Protein (%)	Nitrogen analyzer dry weight basis	≥ 96
Purity (%)	Agarose zone electrophoresis	≥ 96
pH	10% solution	5.0 – 5.6
Moisture (%)	Loss on drying	max. 5.0
Ash (%)	Residue on ignition	≤ 3.0
Storage (°C)		+2 to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.