

<b>Cacodylsäure-Na-Salz·3H<sub>2</sub>O</b> reinst	<b>Kat.Nr./Cat.No. :</b>	<b>15540</b>
Cacodylic acid·Na-salt·3H <sub>2</sub> O research grade		

Parameter parameter	Methode method	Spezifikation specification
<b>Summenformel</b> formula		C <sub>2</sub> H <sub>6</sub> AsO <sub>2</sub> ·Na·3H <sub>2</sub> O
<b>Molekulargewicht</b> molecular weight		214.0
<b>Aussehen</b> appearance		weiße Kristalle white crystals
<b>Gehalt (%)</b> assay	Titration	min. 98.0
<b>Löslichkeit</b> solubility	20% in H <sub>2</sub> O	klar und farblos clear and colorless
<b>Chlorid (ppm)</b> chloride	Cl <sup>-</sup>	max. 200
<b>Gesamtarsen (%)</b> total arsenic	As	min. 34.3
<b>Anorganisches Arsen (%)</b> inorganic arsenic		max. 0.0500
<b>Lagerung (°C)</b> storage		+15 bis / to +30

<p><b>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.</b> We do not guarantee that the product can be used for a special application.</p>
--