
Hefe-Extrakt SERVABACTER gepulvert Yeast extract SERVABACTER powder	Kat.Nr./Cat.No. :	24540
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Aussehen appearance		hellbeiges Pulver light beige powder
Löslichkeit solubility	2% in H ₂ O	klare gelbe Lösung clear yellow solution
pH	2% in H ₂ O	7.0 ± 0.5
Trocknungsverlust (%) loss on drying		max. 5.0
Gesamt-Stickstoff (%) total nitrogen		8.0 – 10.0
Amino-Stickstoff (%) amino nitrogen		2.3 – 4.0
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.