

---

<b>Hefe-Extrakt SERVABACTER®</b> gepulvert Yeast extract SERVABACTER® powder	<b>Kat.Nr./Cat.No. :</b>	<b>24540</b>
---	--------------------------	--------------

---

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification
<b>Aussehen</b> appearance		hellbeiges Pulver light beige powder
<b>Löslichkeit</b> solubility	2 % in H <sub>2</sub> O	klare, gelbe Lösung clear, yellow solution
<b>pH</b>	2 % in H <sub>2</sub> O	7.0 ± 0.5
<b>Trocknungsverlust (%)</b> loss on drying		max. 5.0
<b>Gesamt-Stickstoff (%)</b> total nitrogen		8.0 – 10.0
<b>Amino-Stickstoff (%)</b> amino nitrogen		2.3 – 4.0
<b>Lagerung (°C)</b> storage		+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
We do not guarantee that the product can be used for a special application.