

<b>Kristallviolett</b> reinst Crystal Violet research grade, C.I. 42555	<b>Kat.Nr./Cat.No. :</b>	<b>27335</b>
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
<b>Summenformel</b> formula		C <sub>25</sub> H <sub>30</sub> ClN <sub>3</sub>
<b>Molekulargewicht</b> molecular weight		408.0
<b>Aussehen</b> appearance		grünes Pulver green powder
<b>λ<sub>max</sub> (nm)</b>	0.0001% in H <sub>2</sub> O	586 - 594
<b>A 1cm/ λ<sub>max</sub></b>	0.0001% in H <sub>2</sub> O	min. 0.17
<b>ε<sub>λ<sub>max</sub></sub></b>	0.0001% in H <sub>2</sub> O	min. 69 360
<b>Wasser (%)</b> water	KF	max. 10.0
<b>TLC</b>	Kieselgel 60 F <sub>254</sub> n-Butanol Eisessig/ glacial acetic acid Wasser/ water	einheitlich R <sub>f</sub> ca. 0.8 one spot 4 2 2
<b>Lagerung (°C)</b> storage		+15 bis / to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
We do not guarantee that the product can be used for a special application.