

Product of the Month

SERVA
Electrophoresis

Get **20%** discount* now!

November

*For all SERVA customers within the EU. Outside Germany and Austria online shop orders only. Can not be combined with other discounts, promotions etc.

Protease Inhibitor Mixes

Effective protection against protein degradation

Optimized for protein isolation from various sources

During preparation of cell extracts, proteases are inevitably released from bacteria, yeast, tissue or cell culture. SERVA Protease Inhibitor Mixes provide optimal protection of proteins against proteolytic degradation.

Features

- Easy handling
- High solubility
- Application-optimized
- One vial yields 1 ml of 100x concentrate



Protease Inhibitor Mix G

- For general applications
- Contains 5 protease inhibitors
- Soluble in water

Protease Inhibitor Mix M

- For extracts from mammalian cells
- Contains 6 protease inhibitors
- DMSO for resuspension included

Protease Inhibitor Mix FY

- For fungi and yeast extracts
- Contains 4 protease inhibitors
- DMSO for resuspension included

Protease Inhibitor Mix HP

- Protection of recombinant His-tagged proteins expressed in *E. coli*
- Contains 4 protease inhibitors
- Soluble in water

Ordering information

Product	Cat. no.	Size	List price	Special price*
Protease Inhibitor Mix G	39101.02	5 x 1 vial	113.00 €	90.40 €
	39101.03	10 x 1 vial	201.00 €	160.80 €
Protease Inhibitor Mix M	39102.02	5 x 1 vial	195.00 €	156.00 €
	39102.03	10 x 1 vial	325.00 €	260.00 €
Protease Inhibitor Mix FY	39104.02	5 x 1 vial	154.50 €	123.60 €
	39104.03	10 x 1 vial	283.00 €	226.40 €
Protease Inhibitor Mix HP	39106.02	5 x 1 vial	113.00 €	90.40 €
	39106.03	10 x 1 vial	201.00 €	160.80 €

* The special prices are valid until November, 30th 2016.

For further information, please contact:

Dr. Judith Koch
Tel.: +49 6221 13840-44
E-Mail: Tech.Service@serva.de

SERVA Electrophoresis GmbH
Carl-Benz-Str. 7
D-69115 Heidelberg · Germany
www.serva.de