
| | | |
|---|-------------------|--------------|
| Trypsin NB from porcine pancreas; EC 3.4.21.4 | Cat. No. : | 37294 |
|---|-------------------|--------------|

| Parameter | Method | Specification |
|------------------------------------|--------------------------------|----------------------|
| Activities | Tryptic activity | • 4 000 FIP-U/g |
| Microbiological examination | Colony Forming Units (Ph.Eur.) | Status |

For customer information:

| | |
|---------------------------|--------------------|
| Origin | Porcine pancreas |
| Molecular weight | 24 kDa (trypsin) |
| Appearance | light brown powder |
| Minimum shelf life | 24 months |
| Storage conditions | +2 to +8 °C |

Unit definition:

Tryptic activity (FIP)

1 FIP-U catalyzes the hydrolysis of 1 µmol N-benzoyl-L-arginine ethyl ester (BAEE) per minute at 25°C, pH 8.0.

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.