

SERVA

Market Square

September – December
2016

SDS PAGE

up to
35%



Western Blotting

up to
30%



up to
45%



DNA/RNA
Electrophoresis

Gel Documentation

up to
20%



Relaunched:
www.serva.de

SDS PAGE

Western Blotting

DNA/RNA Electrophoresis

Gel Documentation

The NEW SERVA Website

Precast Gels and Reagents for SDS PAGE

PRiME™ bundles for mini gel SDS PAGE


- Precast gels (cassette 10 cm x 10 cm x 0.7 cm) for high resolution SDS PAGE
- Shelf life up to 12 months for forward-looking stock planning
- Precast gels, buffers and vertical electrophoresis tank – buy packages and save money

- 25%

Cat. no.	Product	Size	Price	Promotion
43296.01	SERVAge™ Package 20 Buffer, consists of: 2x 10 SERVAge™ TG PRiME™ gels of your own choice*, 2 L Laemmli running buffer (10x) plus 20 ml Laemmli sample buffer (2x)	1 kit	309,00 €	231,75 €
43297.01	SERVAge™ Package 20 Chamber, consists of: 2x 10 SERVAge™ TG PRiME™ gels of your own choice*, BlueVertical PRiME™ vertical chamber	1 kit	990,00 €	742,50 €
43283.01	SERVAge™ Package 20 Total, consists of: 2x 10 SERVAge™ TG PRiME™ of your own choice*, 2 L Laemmli running buffer (10x) plus 20 ml Laemmli sample buffer (2x), BlueVertical PRiME™ vertical chamber	1 kit	1.059,00 €	794,25 €



*precast gel type	15 wells	12 wells	10 wells	2D well
SERVAge™ TG PRiME™ 8 %	43284.01	43260.01	43261.01	-
SERVAge™ TG PRiME™ 10 %	43285.01	43263.01	43264.01	-
SERVAge™ TG PRiME™ 12 %	43286.01	43266.01	43267.01	43268.01
SERVAge™ TG PRiME™ 14 %	43287.01	43269.01	43270.01	43271.01
SERVAge™ TG PRiME™ 4 - 12 %	43288.01	43273.01	43274.01	-
SERVAge™ TG PRiME™ 4 - 20%	43289.01	43276.01	43277.01	-
SERVAge™ TG PRiME™ 8 - 16 %	43290.01	43279.01	43280.01	43281.01




Protein standards

- Protein standards ready-to-use, unstained or prestained

- 30%

Cat. no.	Product	Size	Price	Promotion
39248.01	SERVA Unstained Protein Standard II	500 µl	85,00 €	59,50 €
39257.01	SERVA TripleColor Protein Standard II	500 µl	139,00 €	97,30 €
39258.01	SERVA TripleColor Protein Standard III	500 µl	139,00 €	97,30 €
39259.01	SERVA Pink Color Protein Standard II	500 µl	139,00 €	97,30 €




Quick Coomassie® Stain

- Stain your protein gel in one step in 15 minutes only
- For one year stable at room temperature

- 20%

Cat. no.	Product	Size	Price	Promotion
35081.01	Quick Coomassie® Stain	1 L	125,00 €	100,00 €



SDS Pellets

- Due to pellet form of SDS weigh without generating harmful dust

- 20%

Cat. no.	Product	Size	Price	Promotion
20765.03	Dodecylsulfate-Na-Salt in Pellets	1 kg	58,00 €	46,40 €



Acrylamide solutions

- Highest quality from the leading manufacturer of precast polyacrylamide gels in Germany
- No contact with toxic acrylamide powder

- 35%

Cat. no.	Product	Size	Price	Promotion
10680.02	Acrylamide/Bis Solution 40 % w/v 29:1	4 x 500 ml	100,00 €	65,00 €
10681.02	Acrylamide/Bis Solution 40 % w/v 37,5:1	4 x 500 ml	100,00 €	65,00 €
10687.02	Acrylamide/Bis Solution 30 % w/v 29:1	4 x 500 ml	93,00 €	60,45 €
10688.02	Acrylamide/Bis Solution 30 % w/v 37,5:1	4 x 500 ml	93,00 €	60,45 €



SERVA

How to order: phone: +49 6221 13840-47, fax: -10, e-mail: order@serva.de or visit the SERVA online shop at www.serva.de. For SERVA General Conditions of Sale and Supply (GCS) please refer to www.serva.de/gcs. Promotion is valid until December 31st, 2016. Freight charges may apply.

Reagents and Devices for Western Blotting

VisiBlot protein standard

- Seven unstained protein bands featuring multiple IgG binding sites for direct antibody detection
- Three pre-stained protein bands for direct proof of blotting efficiency

- 30%

Cat. no.	Product	Size	Price	Promotion
39260.01	VisiBlot Standard I	500 µl	225,00 €	157,50 €



Xpress Blotting Kit

- Fast and efficient transfer in 15 minutes only
- Ready-to-use – no preparing of buffer or cutting of blotting papers or membranes
- Blotting Fleece instead of several layers of blotting paper for undisturbed transfer

- 30%

Cat. no.	Product	Size	Price	Promotion
42662.01	Xpress Blotting Kit (80 mm x 85 mm)	1 kit	57,00 €	39,90 €
42663.01	Xpress NC Blotting Kit (80 mm x 85 mm)	1 kit	131,00 €	91,70 €
42664.01	Xpress PVDF Blotting Kit (80 mm x 85 mm)	1 kit	179,00 €	125,30 €




SERVALight HRP chemiluminescence kits

- Use SERVALight Vega for all standard ECL applications without the need of changing your existing Western blot protocol
- SERVALight Eos and SERVALight Helios for high demands on signal duration or signal intensity
- Family of chemiluminescence substrate kits for versatile needs

- 20%

Cat. no.	Product	Feature	Size	Price	Promotion
42588.01			50 ml	30,00 €	24,00 €
42588.02	SERVALight Vega	For standard applications	250 ml	114,00 €	91,20 €
42588.03			500 ml	178,00 €	142,20 €
42585.01			50 ml	72,00 €	57,60 €
42585.02	SERVALight Eos	Particularly long standing signal	250 ml	181,00 €	144,80 €
42585.03			500 ml	345,00 €	276,00 €
42587.01			20 ml	95,00 €	76,00 €
42587.02	SERVALight Helios	Particularly high signal strength	100 ml	407,00 €	325,60 €
42587.03			200 ml	649,00 €	519,20 €



BlueBlot Semi-Dry Blotter

- Platinum covered steel net cathode and anode electrodes for homogeneous protein transfer
- Spring mounted anode for flexible applications
- Exchangeable electrode sets (11 cm x 11 cm and 17 cm x 17 cm)
- Stable acryl housing, device and electrodes are very easy to clean

- 20%

Cat. no.	Product	Blotting area	Size	Price	Promotion
BB-SD11	BlueBlot Semi-Dry Blotter SD 11	11 x 11 cm	1 piece	1.685,00 €	1.348,00 €
BB-SD17	BlueBlot Semi-Dry Blotter SD 17	17 x 17 cm	1 piece	1.955,00 €	1.564,00 €
BB-SD26	BlueBlot Semi-Dry Blotter SD 26	24 x 26 cm	1 piece	2.775,00 €	2.220,00 €



GeneGnome

- Optimized and only for imaging of chemiluminescence blots
- Highly sensitive CCD camera (cooled, 4 MP) for best results

- 20%

Cat. no.	Product	Size	Price	Promotion
GGNOME-XRQ-NPC	GeneGnome XRQ (4 MP, 16 bit system without PC, without monitor)	1 system	9.950,00 €	7.960,00 €



SERVA

How to order: phone: +49 6221 13840-47, fax: -10, e-mail: order@serva.de or visit the SERVA online shop at www.serva.de. For SERVA General Conditions of Sale and Supply (GCSS) please refer to www.serva.de/gcss. Promotion is valid until December 31st, 2016. Freight charges may apply.


Reagents and Devices for DNA/RNA Electrophoresis

Agarose SERVA Wide Range

- For agarose gels with a broad separation range from 250 bp up to 23 kb
- For analytical and preparative electrophoresis of nucleic acids
- Excellent solubility

- 45%

Cat. no.	Product	Size	Price	Promotion
11406.02	Agarose SERVA Wide Range	500 g	175,00 €	96,25 €
11406.03		1 kg	335,00 €	184,25 €




SERVA FastLoad DNA standards

- Ready-to-use DNA ladders (50 bp, 100 bp, 1kb), 5 µl/well
- Includes mass/DNA band for easy quantification analysis
- Stable for 12 months at 4 °C, long term storage at -20 °C

- 30%

Cat. no.	Product	Size	Price	Promotion
39315.01	SERVA FastLoad 50 bp DNA Ladder	500 µl	54,00 €	37,80 €
39316.01	SERVA FastLoad 100 bp DNA Ladder	500 µl	54,00 €	37,80 €
39317.01	SERVA FastLoad 1 kb DNA Ladder	500 µl	54,00 €	37,80 €




SERVA DNA Stain Clear G

- Non-carcinogenic alternative to EtBr
- Excitation maximum at 490 nm, additional maxima at 270 nm and 290 nm
- Emission maximum at 530 nm

- 20%

Cat. no.	Product	Size	Price	Promotion
39804.01	SERVA DNA Stain Clear G	1 ml	82,00 €	65,60 €
39804.02		5 x 1 ml	360,00 €	288,00 €



BlueMarine™ submarine chamber

- Available in three formats for a wide application range
- BlueMarine™ HTS for fast analysis of PCR products
- Longlasting, stable gel tanks made from acryl
- Includes gel trays and combs

- 20%

Cat. no.	Product	Gel format	Samples max.	Size	Price	Promotion
BM 100	BlueMarine™ 100	10 cm x 10 cm	28	1 piece	485,00 €	388,00 €
BM 200	BlueMarine™ 200	15 cm x 15 cm und 15 cm x 20 cm	124	1 piece	725,00 €	580,00 €
BM HTS	BlueMarine™ HTS	17.5 cm x 19.2 cm	102	1 piece	1.095,00 €	876,00 €



MP 300 V power supply

- Delivers up to 300 V, 700 mA, 150 W
- Set volt or ampere constant
- 4 x 2 outlets, stackable

- 15%

Cat. no.	Product	Size	Price	Promotion
MP-300V	MP 300 V Power Supply	1 piece	535,00 €	454,75 €



How to order: phone: +49 6221 13840-47, fax: -10, e-mail: order@serva.de or visit the SERVA online shop at www.serva.de. For SERVA General Conditions of Sale and Supply (GCSS) please refer to www.serva.de/gcss. Promotion is valid until December 31st, 2016. Freight charges may apply.

Devices for Gel Documentation

Digital Imaging and Analysis System

- For UV transillumination, white light transillumination and white light top light illumination
- Includes 1D analysis software, not included with the basic version
- UV table, white light table and white light top light are available separately

Cat. no.	Product	Size	Price	Promotion
DIAS-III-B	Digital Imaging and Analysis System III basic, consists of cabinet and camera	1 system	3.495,00 €	2.796,00 €
DIAS-III-L	Digital Imaging and Analysis System III, consists of cabinet and camera and 1D software LabImage L-340	1 system	4.650,00 €	3.720,00 €

- 20%



BlueCube 300

- Equipped with a CCD camera and a two filter system
- UV table (310 nm, 18 cm x 14 cm) is integrated within the drawer
- Includes 1D analysis software

Cat. no.	Product	Size	Price	Promotion
BC-300	SERVA BlueCube 300	1 system	2.895,00 €	2.316,00 €
BC-300L	SERVA BlueCube 300L (incl. laptop)	1 system	3.495,00 €	2.796,00 €

- 20%



SERVA BlueLight Table

- Blue light LED transilluminator for gels up to 20 cm x 12 cm
- To detect DNA in agarose and acrylamide gels stained with EtBr, SERVA DNA Stain Clear G etc.
- No risk of UV damage for high quality work experience

Cat. no.	Product	Size	Price	Promotion
BL-T	SERVA BlueLight Table	1 unit	545,00 €	436,00 €

- 20%



BIO-5000 Plus VIS Gel Scanner

- White light LED scanner for high resolution gel and image documentation
- Transmission and reflective mode
- Short scan times, easy handling

Cat. no.	Product	Size	Price	Promotion
BIO-5000P	BIO-5000 Plus VIS Gel Scanner	1 unit	3.250,00 €	2.762,50 €

- 15%



BIO-1000F Fluorescence Gel Scanner

- Blue light LED scanner unit with filter glass, scan area up to 12.7 cm x 17.8 cm
- Gel scanner for fluorescence stains like SERVA DNA Stain Clear G
- Detects DNA down to 0.04 ng per band

Cat. no.	Product	Size	Price	Promotion
BIO-1000F	BIO-1000F Fluorescence Gel Scanner	1 unit	3.250,00 €	2.762,50 €

- 15%



SERVA






How to order: phone: +49 6221 13840-47, fax: -10, e-mail: order@serva.de or visit the SERVA online shop at www.serva.de. For SERVA General Conditions of Sale and Supply (GCS) please refer to www.serva.de/gcs. Promotion is valid until December 31st, 2016. Freight charges may apply.

The NEW SERVA Website – Online since September 2016



At SERVA surfing scientists –
 → Up-to-date information and best prices*!




-  **NEW:** Choose between English or German language
-  **NEW:** Easy to navigate
-  **NEW:** Online shop including order history for easy re-ordering
-  **NEW:** For PC, tablet and Smartphone
-  **Always updated:** Find here the latest information about our products, webinars, workshops, seminars and exhibitions.

***Found!**
Trypsin MS approved

Trypsin MS approved – a new pack size at an even more attractive pricing

- Each lot is MS qc'ed
- High purity, high stability
- Modified by reductive methylation
- Isolated from porcine pancreas

Cat. no.	Product	Size	Price
37286.01	Trypsin MS approved	1x 100 µg	54,00 €



Request a free sample! Contact SERVA at info@serva.de

SERVA

How to order: phone: +49 6221 13840-47, fax: -10, e-mail: order@serva.de or visit the SERVA online shop at www.serva.de. For SERVA General Conditions of Sale and Supply (GCSS) please refer to www.serva.de/gcss. Promotion is valid until December 31st, 2016. Freight charges may apply.