

Phosphatase Inhibitor Mix I, powder	Cat.No. :	39050
	Lot No.:	P120225

Components	Target-Phosphatase
Imidazol	Alkaline Phosphatase
Sodium fluoride	
Sourum muoride	Acid Phosphatase Protein-Phosphatase-2A Protein-Phosphatase-2B Protein-Phosphatase-2C
Sodium molybdate	Acid Phosphatase Phosphoprotein Phosphatase
Sodium ortho vanadate	Protein-Tyrosine-Phosphatase Alkaline Phosphatase
Sodium tartrate	Acid Phosphatase
Storage (°C)	+15 to +30
Minimum shelf life	05/2014
Reconstitution	Dissolve the content of the vial (156 mg) in 1 ml distilled water to obtain a 100x stock solution (suitable for treatment of 100 ml cell extract). After reconstitution, store at $+2^{\circ}$ C to $+8^{\circ}$ C (stable for 2 years).

We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.

## SERVA Electrophoresis GmbH Quality Control

**Printing date**: 23.05.2012

## Dipl.-Ing. (FH) Bernhard Göckel

Christian Monsler

This report has been computer generated and does not contain a signature.

SERVA Electrophoresis GmbH • D-69115 Heidelberg • Carl-Benz-Str. 7 • Tel.: +49 6221 138 40-0 • Fax +49 6221 138 40-10 Geschäftsführer: Dr. Barbara Müller, Dr. Manfred Kurfürst • HRB 336136 AG Mannheim • E-Mail: info.serva.de • http://www.serva.de COMMERZBANK BLZ: 672 400 39 Kto: 190 38 30 00 • IBAN: DE90 6724 0039 0190 3830 00 • BIC: COBADEFFXXX • Ust.ID-Nr.: DE 812 517 285