## CERTIFICATE OF ANALYSIS



| Protease-Inhibitor-Mix FY | Cat.No. : | 39104   |
|---------------------------|-----------|---------|
|                           | Lot No.:  | P130215 |

Cat. number **Components Target Protease** 12745 **AEBSF-HCl** Serine Proteases 21100 E-64 Cysteine Proteases 52682 **Pepstatin A** Aspartic Proteases 32010 1,10-Phenanthroline Metalloproteases Storage (°C) -15 to -25 Minimum shelf life 06/2016

Reconstitution

Dissolve the content of the vial (62.6 mg) in 1 ml DMSO to obtain a 100x stock solution (suitable for the treatment of 100 ml tissue extract). After reconstitution, aliquot and freeze at -20 °C (stable for at least 4 weeks).

We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH Quality Control

**Printing date**: 10.06.2013

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Steinmann