

Protease-Inhibitor-Mix FY	Cat.No. :	39104
	Lot No. :	P150396

Cat. number	Components	Target Protease	
12745	AEBSF-HCl	Serine Proteases	
21100	E-64	Cysteine Proteases	
52682	Pepstatin A	Aspartic Proteases	
32010	Pepstatin A 1,10-Phenanthroline	Metalloproteases	
	×10.		
Storage (°C)	90	-15 to -25	
Minimum shelf life	CU	20.10.2020	
Winning Siten file		20.10.2020	
	IA		
Reconstitution	Dissolve the content of the vial (6	52.6 mg) in 1 ml DMSO to obtain a	
0	100x stock solution (suitable for the treatment of 100 ml tissue extract).		
	After reconstitution, aliquot and freeze at -20 °C (stable for at least 4		
	weeks).		

We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH Quality Control

Printing date: 04.01.2019

Christian Monsler

Patricia Milford

This report has been computer generated and does not contain a signature.