

Protease-Inhibitor-Mix HP	Cat.No. :	39106
	Lot No.:	P130489

Cat. number	Components	Target Protease
12745	AEBSF-HCl	Serine Proteases
13718	Aprotinin from bovine lung	Serine Proteases
21100	E-64	Cysteine Proteases
51867	Leupeptin	Cysteine Proteases and Trypsin-like Serine Proteases
Storage (°C)	C ^{YO}	-15 to -25
Minimum shelf life		12/ 2016
Reconstitution	Dissolve the content of the vial (12.7 m a 100x stock solution (suitable for the	

a 100x stock solution (suitable for the treatment of 100ml tissue extract). After reconstitution, aliquot and freeze at -20 °C (stable for 4 weeks).

We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.

SERVA Electrophoresis GmbH Quality Control **Printing date**: 04.12.2013

Dipl.-Ing. (FH) Bernhard Göckel

Patricia Milford

This report has been computer generated and does not contain a signature.

SERVA Electrophoresis GmbH • D-69115 Heidelberg • Carl-Benz-Str. 7 • Tel.: +49 6221 138 40-0 • Fax +49 6221 138 40-10 Geschäftsführer: Dr. Barbara Müller, Dr. Manfred Kurfürst • HRB 336136 AG Mannheim • E-Mail: info.serva.de • http://www.serva.de COMMERZBANK BLZ: 672 400 39 Kto: 190 38 30 00 • IBAN: DE90 6724 0039 0190 3830 00 • BIC: COBADEFFXXX • Ust.ID-Nr.: DE 812 517 285