

Protease-Inhibitor-Mix HP	Cat.No. :	39106
	Lot No.:	P140548

Cat. number	Components	Target Protease	
		$\sim$	
12745	AEBSF-HCI	Serine Proteases	
13718	Aprotinin from bovine lung	Serine Proteases	
21100	E-64	Cysteine Proteases	
51867	Leupeptin	Cysteine Proteases and Trypsin-like Serine Proteases	
Storage (°C)	ctrop	-15 to -25	
Minimum shelf life		12/ 2017	
Reconstitution	Dissolve the content of the vial (12.7 r		
a 100x stock solution (suitable for the treatment of 100ml tissue extract).			

After reconstitution, aliquot and freeze at -20 °C (stable for 4 weeks).

We do not guarantee that the product can be used for a special application. This document does not release you from performing the standard control upon receipt of incoming goods.

**SERVA Electrophoresis GmbH Quality Control** 

**Printing date**: 09.12.2014

Dipl.-Ing. (FH) Bernhard Göckel

Daniela Lux-Helmstetter

This report has been computer generated and does not contain a signature.

SERVA Electrophoresis GmbH • D-69115 Heidelberg • Carl-Benz-Str. 7 • Tel.: +49 6221 138 40-0 • Fax +49 6221 138 40-10 Geschäftsführer: Dr. Barbara Müller, Dr. Manfred Kurfürst • HRB 336136 AG Mannheim • E-Mail: info.serva.de • http://www.serva.de COMMERZBANK BLZ: 672 400 39 Kto: 190 38 30 00 • IBAN: DE90 6724 0039 0190 3830 00 • BIC: COBADEFFXXX • Ust.ID-Nr.: DE 812 517 285