

Acarbose from Actinoplanes sp.	Kat.Nr./Cat.No. :	10017
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₂₅ H ₄₃ O ₁₈ N
Molekulargewicht molecular weight		645.6 Da
Aussehen appearance		weißes bis gelbliches Pulver white to yellowish powder
Gehalt (%) assay	HPLC	95.0 – 101.0
Verunreinigungen (%) sum of unspecified impurities		max. 1.0
Mikrobiologische Reinheit (CFU/g) microbial purity total	total aerobic microbial count (TAMC)	max. 1000
	total combined yeast/ mould count (TYMC)	max. 100
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.