

Acrylamid/ Bis-Lösung, 37.5 : 1 (40%w/v) Acrylamide/ Bis Solution, 37.5 : 1 (40%w/v)	Kat.Nr./Cat.No. :	10681
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Beschreibung description	Lösung enthält 40% (w/v) Acrylamid-Bisacrylamid (37.5:1) in deionisiertem Wasser. Solution contains 40% (w/v) acrylamide-bisacrylamide (37.5:1) in deionized water.	
Aussehen appearance		klar und farblos clear and colorless
Gehalt (%) assay	aus Brechungsindex from refractive index	40.0 ± 0.5
pH	5% in H ₂ O	6.0 – 8.0
A 1cm/ 290 nm	5% in H ₂ O	max. 0.7
Acrylsäure (%) acrylic acid	Titration	max. 0.03
Leitfähigkeit (µS/cm) conductivity		max. 100
Polymerisationszeit (min) polymerizing time	15% Gel / 25°C	max. 30
Lagerung (°C) storage		+2 bis/ to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.