

Adenosin-5'-diphosphat·Na₂-salz (ADP) reinst

Kat.Nr./Cat.No. : 10800

Adenosine-5'-diphosphate·Na₂-salt (ADP) research grade

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₁₀ H ₁₃ N ₅ O ₁₀ P ₂ ·Na ₂
Molekulargewicht molecular weight		471.2
Aussehen appearance		weißes Pulver white powder
ADP·Na₂ (%)	enzym. calc. on dry basis	min. 96.0
ADP (%)	enzym.	min. 82.0
AMP (%)	enzym.	max. 3.0
ATP (%)	enzym.	max. 1.0
Wasser (%) water	KF	max. 7.0
Lagerung (°C) storage		-15 bis/ to -25

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.