

Adenosin-5'-phosphat-Na₂-Salz (AMP) p.A.	Kat.Nr./Cat.No. :	10883
Adenosine-5'-phosphate-Na ₂ -salt (AMP), analytical grade		

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula	AMP•Na ₂ AMP•Na ₂ •7H ₂ O	C ₁₀ H ₁₂ N ₅ O ₇ P•Na ₂ C ₁₀ H ₁₂ N ₅ O ₇ P•Na ₂ •7H ₂ O
Molekulargewicht molecular weight		391.2 (wasserfrei/ anhydrous) 517.2
Aussehen appearance		weißes Pulver white powder
Gehalt (%) assay	HPLC	min. 98.0
Wasser (%) water	KF	max. 26.0
Schwermetalle (ppm) heavy metals	als/ as Pb	max. 10
Ursprung origin		nicht tierisch non animal
Lagerung (°C) storage		+2 bis / to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.