
Adenosin-5'-triphosphat·Na₂-Salz (ATP) krist. reinst	Kat.Nr./Cat.No. :	10920
Adenosine-5'-triphosphate·Na ₂ -salt (ATP) cryst. research grade		

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₁₀ H ₁₄ N ₅ O ₁₃ P ₃ ·Na ₂
Molekulargewicht molecular weight		551.1
Aussehen appearance		weißes Pulver white powder
Gehalt (%) assay	HPLC	min. 98.0
Wasser (%) water	KF	max. 10.0
Lagerung (°C) storage		-15 bis/ to -25

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.