
Ethanol unvergällt 96% p.A. Ethanol undenaturated 96% analytical grade	Kat.Nr./Cat.No. :	11094
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₂ H ₆ O
Molekulargewicht molecular weight		46.07
Aussehen appearance		klare farblose Flüssigkeit clear colorless liquid
Wasser (%) water	KF	4.0 – 6.0
Brechungsindex refractive index	n _D /20°C	1.3625 – 1.3635
Identität identity	IR	entspricht corresponds
Lagerung (°C) storage		+15 bis/ to +30

SERVA Electrophoresis GmbH

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.