

Ethylendiamin-tetraessigsäure· K₃-salz p.A.	Kat.Nr./Cat.No. :	11277
Ethylenediamine tetraacetic acid· K ₃ -salt analytical grade		

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₁₀ H ₁₃ N ₂ O ₈ ·K ₃ ·2H ₂ O
Molekulargewicht molecular weight		442.55
Aussehen appearance		weißes krist. Pulver white cryst. powder
Gehalt (%) assay	Titration	min. 99.0
A 1cm/ 260 nm	5% in H ₂ O	max. 1.7
A 1cm/ 280 nm	5% in H ₂ O	max. 0.5
pH	5% in H ₂ O	7.0 – 9.0
Wasser (%) water	KF	7.5 – 9.0
Lagerung (°C) storage		+15 bis / to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.