

Agar Agar SERVA Kobe 1 Pulver reinst Agar Agar SERVA Kobe 1 powder research grade	Kat.Nr./Cat.No. : 11392
---	---------------------------------------

Parameter parameter	Methode method	Spezifikation specification
Aussehen appearance		weißes bis gelbliches Pulver white to yellowish powder
Gelstärke (g/cm²) gel strength	1.5% Gel 1.5% gel	min. 800
Geliertemperatur (°C) point of solidification	1.5% Lösung 1.5% solution	35 – 43
A 1cm/430nm A 1cm/525nm	1.5% in H ₂ O (60°C) 1.5% in H ₂ O (60°C)	max. 0.5 max. 0.45
Trocknungsverlust (%) loss on drying	4h, 105°C	max. 15.0
Asche (%) ash		max. 1.5
pH	1.5% in H ₂ O (60°C)	6.0 – 8.0
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.