

Agar Agar SERVA Kobe 1 in Streifen, reinst Agar Agar SERVA Kobe 1 in stripes, research grade	Kat.Nr./Cat.No. : 11395
--	---------------------------------------

Parameter parameter	Methode method	Spezifikation specification
Aussehen appearance		hellelfenbeinfarbene Streifen light ivory-colored strips
Gelstärke (g/cm²) gel strength	1.5 % Gel 1.5 % gel	> 400
Gelertemperatur(°C) point of solidification	1.5 % Lösung 1.5 % solution	30 – 40
A 1 cm/ 430 nm	1.5 % in H ₂ O (60 °C)	max. 0.5
A 1 cm/ 525 nm	1.5 % in H ₂ O (60 °C)	max. 0.45
Trocknungsverlust (%) loss on drying	4 h, 105 °C	max. 25.0
Asche (%) ash		max. 6.5
pH	1.5 % in H ₂ O (60 °C)	5.0 – 8.0
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.