

Agarose SERVA hohe EEO Agarose SERVA high EEO	Kat.Nr./Cat.No. : 11397
---	---------------------------------------

Parameter parameter	Methode method	Spezifikation specification
Aussehen appearance		reinweißes Pulver white powder
Gelstärke (g/cm²) gel strength	1.5% Gel 1.5% gel	min. 1300
Geliertemperatur (°C) point of solidification	1.5% Lösung 1.5% solution	34 – 39
A 1cm/430nm A 1cm/525nm	1.5% in H ₂ O (60°C) 1.5% in H ₂ O (60°C)	max. 0.1 max. 0.05
Trocknungsverlust (%) loss on drying	4h, 105°C	max. 10.0
Asche (%) ash		max. 1.0
Sulfat (%) sulphate		max. 0.6
Elektroendosmose (mr) electroendosmosis		0.23 – 0.27
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.