

Agarose SERVA neutral für IEF Agarose SERVA neutral for IEF.	Kat.Nr./Cat.No. :	11402
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Aussehen appearance		weißes Pulver white powder
Gelstärke (g/cm²) gel strength	1.5% Gel 1.5% gel	> 800
Geliertemperatur (°C) point of solidification	1.5% Lösung 1.5% solution	38 – 43
A 1cm/430nm A 1cm/525nm	1.5% in H ₂ O (60°C) 1.5% in H ₂ O (60°C)	max. 0.15 max. 0.02
Trocknungsverlust (%) loss on drying	4h, 105°C	max. 5.0
Asche (%) ash		max. 0.7
Elektroendosmose (mr) electroendosmosis		± 0.0
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.