

Agarose SERVA für DNA-Elektrophorese reinst Agarose SERVA for DNA-electrophoresis research grade	Kat.Nr./Cat.No. : 11404
--	---------------------------------------

Parameter parameter	Methode method	Spezifikation specification
Aussehen appearance		weißes Pulver white powder
Gelstärke (g/cm²) gel strength	1.5% Gel	min. 1700
Gelertemperatur (°C) point of solidification	1.5% Lösung 1.5% solution	34 – 39
A 1cm/ 430nm	1.5% in H ₂ O, 60°C	max. 0.05
A 1cm/ 525nm	1.5% in H ₂ O, 60°C	max. 0.02
Trocknungsverlust (%) loss on drying	4h, 105°C	max. 8.5
Asche (%) ash		max. 0.7
Elektroendosmose (-Mr) electroendosmosis		max. 0.13
Sulfat (%) sulphate		max. 0.14
DNA-Elektrophorese DNA-electrophoresis		entspricht corresponds
Ursprung origin		aus reinem Seetang from raw seaweeds
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.