

Agarose SERVA Tabletten, 0.5 g/Tablette molecular biology grade Agarose SERVA tablets, 0.5 g/tablet	Kat.Nr./Cat.No. :	11405
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Aussehen appearance		weiße Tabletten white tablets
Gelstärke (g/cm²) gel strength	1.5% Gel	> 1000
Geliertemperatur (°C) point of solidification	1.5% Lösung 1.5% solution	34 – 39
A 1cm/430nm	1.5% in H ₂ O, 60°C	max. 0.05
A 1cm/525nm	1.5% in H ₂ O, 60°C	max. 0.02
Trocknungsverlust (%) loss on drying	4h, 105°C	max. 8.0
Asche (%) ash		max. 0.7
Elektroendosmose (-Mr) electroendosmosis		max. 0.13
Sulfat (%) sulphate		max. 0.1
DNA-Elektrophorese DNA-electrophoresis		entspricht corresponds
Ursprung origin		aus reinem Seetang from raw seaweeds
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.