

Albumin aus Hühnerei lyophil. salzfrei Albumin egg (ovalbumin) lyophil. salt-free	Kat.Nr./Cat.No. :	11840
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Molekulargewicht molecular weight		ca. 45 000
Aussehen appearance		weißes Lyophilisat white lyophilisate
Löslichkeit solubility	1% in PBS	klar bis leicht trübe Lösung clear to slightly turbid solution
Wasser (%) water	K F	max. 4.0
Gehalt (%) (Ovalbumin) assay	SDS-PAGE 15%	min. 80.0
Conalbumin (%)	SDS-PAGE 15%	max. 7.0
Ovomucoid (%)	SDS-PAGE 15%	max. 1.0
Lyozym (%) lyozyme	SDS-PAGE 15%	nicht nachweisbar not detectable
Salmonellen salmonella		nicht nachweisbar not detectable
Lagerung (°C) storage		+2 bis/ to +8

<p>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert. We do not guarantee that the product can be used for a special application.</p>
--