

**Albumin bovine**  
cryst. lyophil.

**Cat.No. : 11920**

Parameter	Method	Specification
<b>Molecular weight</b>		ca. 67 000
<b>Protein (%)</b> (protein factor 6.22)	Nitrogen Analyzer dry weight basis	98.0 – 100.0
<b>Purity (%)</b> (albumin, 7% solution)	Cellulose Acetate Electrophoresis	99.0 – 100.0
<b>pH</b>	7% solution	5.0 – 5.4
<b>Sodium (mg/g)</b>	Flame Photometer	0 – 5.0
<b>Chloride (mg/g)</b>	Coulometric Titrator	0 – 3.0
<b>Moisture (%)</b>	Karl Fischer	0 – 5.0
<b>Storage (°C)</b>		+2 bis/ to +8

SERVA Electrophoresis GmbH

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
We do not guarantee that the product can be used for a special application.