

Albumin bovine Fraction V, pH 5.2
standard grade, lyophil.

Cat.No.: 11922

Parameter	Method	Specification
Molecular weight		ca. 67 000
Appearance		white lyophilisate
Protein (%) (protein factor 6.22)	Nitrogen Analyzer dry weight basis	96 – 100
Purity (%) (albumin, 7% solution)	Cellulose Acetate Electrophoresis	98 – 100
pH	1% solution	5.2 – 5.6
Sodium (mg/g)	Flame Photometer	0 – 5.0
Chloride (mg/g)	Coulometric Titrator	0 – 3.0
Moisture (%)	KF	0 – 5.0
Storage (°C)		+2 bis +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.