

Albumin bovine Fraction V, pH 7.0
standard grade, lyophil.

Cat.No. : 11930

Parameter	Method	Specification
Molecular weight		ca. 67 000
Protein (%) (protein factor 6.22)	Nitrogen Analyzer dry weight basis	96 – 100
Purity (%) (albumin, 7% solution)	Cellulose Acetate Electrophoresis	98 – 100
pH	7% solution	6.8 – 7.2
Sodium (mg/g)	Flame Photometer	0 – 10.0
Chloride (mg/g)	Coulometric Titrator	0 – 3.0
Iron (µg/g)	Atomic Absorption	0 – 15.0
IgG	Radialimmunodiffusion	not detectable
Moisture (%)	KF	0 – 5.0
Storage (°C)		+2 to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.