

**Albumin bovine Cohn Fraction V pH 7**  
standard grade, lyophil.

**Cat.No. : 11931**

Parameter	Method	Specification
<b>Molecular weight</b>		ca. 67 000
<b>Protein (%)</b> (protein factor 6.25)	Nitrogen Analyzer dry weight basis	≥ 96
<b>Solubility</b>	4% solution, • A550	≤ 0.1
<b>Color</b>	4% solution, • A404	≤ 0.4
<b>Water (%)</b>	Loss on drying	≤ 6.0
<b>Sulphated ash (%)</b>		≤ 3.0
<b>pH</b>	1% solution in 0.15 M NaCl	6.8 – 7.2
<b>Heavy metals (ppm)</b>	AAS	≤ 50
<b>Albumin (%)</b>	Electrophoresis	≥ 96.0
<b>Storage (°C)</b>		+2 to +8

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
We do not guarantee that the product can be used for a special application.