

Albumin bovine Fraction V
Fatty Acid-free, lyophil.

Cat.No.: 11932

Parameter	Method	Specification
Molecular weight		ca. 67 000
Appearance		off-white lyophilisate
Protein (%) (protein factor 6.22)	Nitrogen Analyzer dry weight basis	96 – 100
Purity (%) (albumin, 7% solution)	Cellulose Acetate Electrophoresis	98 – 100
pH	7 % solution	6.8 – 7.2
Sodium (mg/g)	Flame Photometer	0 – 10.0
Chloride (mg/g)	Coulometric Titrator	0 – 6.0
Total lipids (mg/g)	Supercritical Fluid Extraction	0 – 2.0
Free fatty acids (mg/g)	Chromatography	0 – 0.2
Esterified fatty acids (mg/g)	Chromatography	0 – 1.0
Moisture (%)	KF	0 – 5.0
Storage (°C)		+2 to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.