

Albumin bovine Fraction V, pH 7.0
microbiological grade, lyophil.

Cat.No.: 11934

| Parameter | Method | Specification |
|---|---------------------------------------|----------------------------|
| Molecular weight | | ca. 67 000 |
| Appearance | | ivory colored lyophilisate |
| Protein (%) (protein factor 6.22) | Nitrogen Analyzer dry weight basis | 96 – 100 |
| Purity (%) (albumin, 7% solution) | Cellulose Acetate Electrophoresis | 98 – 100 |
| pH | 7% solution | 6.8 - 7.2 |
| Moisture (%) | Karl Fischer | 0 – 5.0 |
| Storage (°C) | | +2 to +8 |

SERVA Electrophoresis GmbH

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.