

Albumin Bovine 30% solution
polymer enhanced

Cat.No.: 11937

Parameter	Method	Specification
Protein (g/dL)	Biuret	29 - 31
pH		7.2 - 7.4
Sodium Chloride (g/dL)	Coulometric Titrator	0.60 - 0.70
IgG	Radialimmunodiffusion	not detectable
Immunohematological Performance	-Titration of Rh system antibody -Hemagglutination slide test with Rh system antibody -Microscopic examination of RBC morphological characteristics -T-Activation of RBCs -Screening against panel of known RBCs at room temperature, 37°C, high protein phase, and indirect antiglobulin test	complies
Preservative (g/dL)	Sodium Azide	≤ 0.1
Storage (°C)		+2 to +8

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.