

<b>Amidoschwarz 10 B</b> Amido Black 10 B	<b>Kat.Nr./Cat.No. :</b>	<b>12310</b>
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
<b>Summenformel</b> formula		$C_{22}H_{14}N_6O_9S_2 \cdot Na_2$
<b>Molekulargewicht</b> molecular weight		616.5
<b>Aussehen</b> appearance		schwarzbraunes Pulver black to brown powder
<b><math>\lambda_{max}</math> (nm)</b>	0.001% in H <sub>2</sub> O	616 - 624
<b>A 1cm/<math>\lambda_{max}</math></b>	0.001% in H <sub>2</sub> O	min. 0.4
<b><math>\epsilon_{\lambda_{max}}</math></b>	0.001% in H <sub>2</sub> O	min. 24 660
<b>Wasser (%)</b> water	KF	max. 15.0
<b>TLC</b>	Kieselgel/silicagel 60 F <sub>254</sub> n-Propanol 4 i-Butanol 2 Ethylacetat 1 Wasser/water 3	1 Hauptkomponente RF: ca. 0.7 Spuren möglich one main spot, traces possible
<b>Lagerung (°C)</b> storage		+15 bis/ to +30

<p><b>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.</b> We do not guarantee that the product can be used for a special application.</p>
--