

| | | |
|---|--------------------------|--------------|
| AEBSF • HCl reinst 4-(2-Aminoethyl)-benzene sulfonyl fluoride-HCl research grade | Kat.Nr./Cat.No. : | 12745 |
|---|--------------------------|--------------|

| Parameter parameter | Methode method | Spezifikation specification |
|---|--------------------------|---|
| Summenformel formula | | C ₈ H ₁₀ FNO ₂ S•HCl |
| Molekulargewicht molecular weight | | 239.69 |
| Aussehen appearance | | weißes Pulver white powder |
| Gehalt (%) assay | HPLC | min. 98.0 |
| Schmelzpunkt (°C) melting point | | 181 – 186 |
| Löslichkeit solubility | 20% in H ₂ O | klar und farblos clear and colorless |
| | 1% in MeOH | klar und farblos clear and colorless |
| Lagerung (°C) storage | | +2 bis/ to +8 |

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.