

# SPEZIFIKATION

## SPECIFICATION

---

<b>Ammoniumformiat</b> reinst Ammonium formate research grade	<b>Kat.Nr./Cat.No. :</b>	<b>13369</b>
--	--------------------------	--------------

---

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification
<b>Summenformel</b> formula		HCOONH <sub>4</sub>
<b>Molekulargewicht</b> molecular weight		63.06
<b>Aussehen</b> appearance		farblose Kristalle colorless crystals
<b>Gehalt (%)</b> assay	Titration	min. 98.0
<b>pH</b>	5% in H <sub>2</sub> O	6.0 – 7.0
<b>Schmelzpunkt (°C)</b> melting point		ca. 116
<b>Schwermetalle (%)</b> heavy metals	Pb	max. 0.002
<b>Chlorid (%) / chloride</b>	Cl	max. 0.005
<b>Sulfat (%) / sulfate</b>	SO <sub>4</sub>	max. 0.01
<b>Eisen (%) / iron</b>	Fe	max. 0.001
<b>Lagerung (°C)</b> storage		