

Ammoniummolybdat p.A. Ammonium molybdate analytical grade	Kat.Nr./Cat.No. :	13370
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$(\text{NH}_4)_6\text{Mo}_7\text{O}_{24}\cdot 4\text{H}_2\text{O}$
Molekulargewicht molecular weight		1235.9
Aussehen appearance		weiße Kristalle white crystals
Gehalt (%) assay	Titration als / as MoO_3 als / as $(\text{NH}_4)_6\text{Mo}_7\text{O}_{24}\cdot 4\text{H}_2\text{O}$	81.0 – 83.0 min. 99.0
Chlorid (%) chloride	Cl^-	max. 0.002
Phosphat (%) phosphate		max. 0.0005
Blei (%) lead	Pb	max. 0.001
Identität identity	IR	entspricht corresponds
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.