

**Anilinblau** pract. **Kat.Nr./Cat.No.:** **13645**  
Aniline Blue pract.

Parameter parameter	Methode method	Spezifikation specification
<b>Summenformel</b> formula		$C_{37}H_{27}N_3Na_2O_9S_3$
<b>Molekulargewicht</b> molecular weight		799.8
<b>Aussehen</b> appearance		rotschwarze Kristalle red black crystals
<b><math>\lambda_{max}</math> (nm)</b>	0.001% in H <sub>2</sub> O	570 - 610
<b>A 1cm/<math>\lambda_{max}</math></b>	0.001% in H <sub>2</sub> O	min. 0.2
<b>Wasser (%)</b> water	KF	max. 10.0
<b>TLC</b>	Kieselgel 60 F <sub>254</sub> n-Butanol Ethanol Eisessig/ glacial acetic acid Wasser/ water	1 Hauptbande 5 2 1 2
<b>Lagerung (°C)</b> storage		+15 bis/ to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
We do not guarantee that the product can be used for a special application.