

9-Anthryldiazomethan reinst **Kat.Nr./Cat.No. : 13682**
9-Anthryldiazomethane research grade

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₁₅ H ₁₀ N ₂
Molekulargewicht molecular weight		218.6
Aussehen appearance		rotorangees krümeliges Pulver reddish orange-colored crumbly powder
λ max (nm)	0.001 % in Aceton 0.001 % in acetone	391 - 399
A 1 cm/ λ max	0.001 % in Aceton 0.001 % in acetone	min. 0.28
Löslichkeit solubility	1 mg/ ml in Aceton 1 mg/ ml in acetone	klare gelbe Lösung clear yellow solution
Lagerung (°C) storage		-15 bis/ to -25

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.