

Antipain reinst research grade	Kat.Nr./Cat.No. :	13685
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$C_{27}H_{44}N_{10}O_6 \cdot HCl \cdot 2H_2O$
Molekulargewicht molecular weight		677.2
Aussehen appearance		schwach gelbes Pulver pale yellow powder
pH	1 % in H ₂ O	5.0 – 6.0
Trocknungsverlust (%) loss on drying		max. 5.0
Aktivität (µg/ml) activity	IC50	max. 0.5
Lagerung (°C) storage		-15 bis/ to -25
Referenzen references	J. Antibiot., 25, 263 - 266 (1972) J. Antibiot., 25, 267 - 270 (1972)	

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.