

Cellulose mikrokristallin ca. 0.02mm reinst Cellulose microcrystalline ca. 0.02mm research grade	Kat.Nr./Cat.No. :	14205
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Aussehen appearance		weißes Pulver white powder
Gehalt (%) assay		97.0 – 102.0
pH	12.5 % in H ₂ O	5.0 – 7.0
Trocknungsverlust (%) loss on drying	3 h, 105°C	max. 6.0
Wasserlösliche Substanzen water soluble substances	mg/ 5g	max. 12.5
Schwermetalle (ppm) heavy metals		max. 10
Test auf Stärke starch test		negativ negative
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.