

Azorubin rein Azorubin pure	Kat.Nr./Cat.No. :	14410
---------------------------------------	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		$C_{20}H_{12}N_2O_7S_2 \cdot Na_2$
Molekulargewicht molecular weight		502.4
Aussehen appearance		rotbraunes Pulver reddish brown powder
Reinheit purity	TLC	einheitlich homogeneous
λ max (nm)	0.001 % in H ₂ O	514 – 522
A 1 cm/ λ max	0.001 % in H ₂ O	min. 0.4
ϵ / λ max	0.001 % in H ₂ O	min. 20 000
Wasser (%) water	KF	max. 15.0
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.