

Bathocuproin disulfonsäure•Na₂-salz p.A. Bathocuproine disulfonic acid•Na ₂ -salz analytical grade	Kat.Nr./Cat.No. : 14470
--	---------------------------------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₂₆ H ₁₈ N ₂ O ₆ S ₂ •Na ₂
Molekulargewicht molecular weight		564.6
Aussehen appearance		beiges Pulver beige powder
Gehalt (%) assay	Titration	min. 98.0
Wasser (%) water	KF	max. 4.0
Identität identity	IR	entspricht corresponds
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.