

6-Benzylaminopurin (N⁶-Benzyladenin) reinst 6-Benzylaminopurine (N ⁶ -Benzyladenine) research grade	Kat.Nr./Cat.No. :	14812
---	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₁₂ H ₁₁ N ₅
Molekulargewicht molecular weight		225.25
Aussehen appearance		weiße Kristalle white crystals
Gehalt (%) assay	HPLC	min. 99.0
λ max (nm)	0.001% in 0.1M HCl	274 – 278
Schmelzpunkt (°C) melting point		230 – 233
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.