

Brenztraubensäure-Na-salz reinst Pyruvic acid-Na-salt research grade	Kat.Nr./Cat.No. :	15220
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₃ H ₃ O ₃ ·Na
Molekulargewicht molecular weight		110.04
Aussehen appearance		fast weißes Pulver almost white powder
Gehalt (%) assay	Titration	min. 99.0
Löslichkeit solubility	1.5g in 25ml H ₂ O	klare, leicht gelbliche Lösung clear, light yellowish solution
Trocknungsverlust (%) loss on drying		max. 0.5
Sulfatasche (%) sulfated ash		60.0 – 70.0
Schwermetalle (ppm) heavy metals	als/ as Pb	max. 10
Freie Säure (%) free acid	Titration als/ as C ₃ H ₄ O ₃	max. 1.0
Arsen (ppm) arsenic	As ₂ O ₃	max. 1
Ursprung origin		nicht tierisch non animal
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.