

<b>5-Brom-4-chlor-3-indolyl-β-D-glucuronid-cyclohexylammoniumsalz</b> 5-Bromo-4-chloro-3-indolyl-β-D-glucuronide cyclohexylammonium salt	<b>Kat.Nr./Cat.No. :</b>	<b>15248</b>
--	--------------------------	--------------

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification
<b>Summenformel</b> formula		C <sub>20</sub> H <sub>26</sub> BrClN <sub>2</sub> O <sub>7</sub>
<b>Molekulargewicht</b> molecular weight		521.8
<b>Aussehen</b> appearance		weißes krist. Pulver white cryst. powder
<b>Reinheit</b> purity	TLC	min. 98
<b>Löslichkeit</b> solubility	2 % in DMF	klare farblose Lösung clear colorless solution
<b>Wasser (%)</b> water	KF	max. 3.5
<b>[α]<sub>D</sub>/20°C</b>	1 % in DMF/ H <sub>2</sub> O 1:1	-87.0° bis/ to -93.0°
<b>Identität</b> identity	IR	entspricht corresponds
<b>Lagerung (°C)</b> storage		-15 bis/ to -25

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.  
We do not guarantee that the product can be used for a special application.