

5-Brom-4-chlor-3-indoxylphosphat-Na₂-salz p.A. 5-Bromo-4-chloro-3-indolyl-phosphate·Na ₂ -salt analytical grade	Kat.Nr./Cat.No. :	15259
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₈ H ₄ BrClNO ₄ P·Na ₂
Molekulargewicht molecular weight		370.43
Aussehen appearance		fast weißes bis schwach gelbliches Pulver off-white to pale yellow powder
Gehalt (%) assay	HPLC	min. 99.0
Wasser (%) water	KF	max. 10.0
Identität identity	IR	entspricht corresponds
Lagerung (°C) storage		-15 bis/ to -25

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.